Ambu†

Ambu[®] SPUR[®] II

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Directions for use

These directions for use may be updated without further notice. Copies of the current version are available from the manufacturer.

For use by CPR-trained personnel only.

The resuscitator must be unpacked and prepared for immediate use before it is placed ready for use in emergency situations:


1. If the resuscitator is packed in a compressed state, unfold the resuscitator by pulling on the patient valve and the inlet valve (if the resuscitator is supplied with a bag reservoir, do not pull this as tearing may occur).
2. Make a brief functional check as described in paragraph 8 of this manual.
3. Fit the face mask. Supply pharyngeal airways and other equipment prescribed for the intended standby location and place all items in the plastic bag supplied with the resuscitator. The integrity of kits issued for storage ready for use should be inspected at the interval established by your local protocol.
4. Single patient use resuscitator, no service required.

CAUTION

Make sure that the resuscitator is never stored in a deformed state other than as folded when delivered by the manufacturer, otherwise permanent distortion of the bag will occur which may reduce the ventilation efficiency. The folding zone is clearly visible on the bag (only Adult and Paediatric versions may be folded).

For illustration purposes, only the Paediatric version is shown in the rest of this Directions for Use, except where specifically stated.

*Also available with pressure limiting valve and manometer port.

1. Warnings / precautions

Failure to observe these precautions may result in inefficient ventilation of the patient or damage to the equipment.

WARNING

- Oil or grease should not be used in close proximity to oxygen equipment - fire may result.
- Do not smoke or use open flames when oxygen is in use - fire may result.
- Never override the pressure-limiting valve (if available) unless medical and professional assessment indicates the necessity. High ventilation pressures may cause lung rupture to certain patients. If the pressure-limiting valve is overridden, a manometer must be used to monitor ventilatory pressures and avoid the possibility of lung rupture. See paragraph 6.1 for further details.
- If an exhalation filter is used: The exhalation filter may create a small amount of inadvertent breathing resistance at 50 l/min. Do not use this device if this additional breathing resistance would be detrimental for the patient. If the exhalation filter becomes obstructed by foreign material or offers excessive resistance to exhalation, discontinue use and discard.

CAUTION

- For use by CPR-trained personnel only. The proper application of a face mask to obtain a tight seal should be trained in particular (paragraph 6.1). Make sure that such personnel are made familiar with the contents of this manual.
- Always perform a functional test of the resuscitator after unpacking and assembly. See section 8 for further information.
- Do not clean or sterilize the Ambu resuscitator, as it is for use on a single patient only. No reprocessing must be done.
- The efficiency of the ventilation should always be checked by watching the movement of the chest and listening for the expiratory flow from the valve.
- Switch immediately to mouth-to-mouth resuscitation if efficient ventilation cannot be obtained with the device.
- Do not use the resuscitator in toxic or hazardous atmospheres.
- Change to injection directly in the tube if unusually high flow resistance is felt through the MediPort.

2. Intended use

The Ambu® SPUR® II resuscitator is intended for pulmonary resuscitation. The range of application for each version is:

- Adult: Adults and children with a body weight more than 30 kg (66 lbs).
- Paediatric: Infants and children with a body weight up to 30 kg (66 lbs), approx. 9 years of age.
- Infant: Neonates and infants with a body weight up to 10 kg (22 lbs).

MediPort allows injection of drug (through the MediPort) directly in the endotracheal tube without disassembling the tube/patient valve. At the same time, the resuscitator bag can be used as normal for continued oxygen supply and assisted ventilation.

3. Specifications

The Ambu SPUR II resuscitator conforms to the following standards: EN ISO 10651-4:2002 and ASTM F 920-93.

The Ambu SPUR II resuscitator conforms to Council Directive 93/42/EEC concerning Medical Devices.

		Infant:	Paediatric:	Adult:
Stroke volume:	One hand	150 ml	450 ml	800 ml
	Two hands	-	-	1100 ml
Resuscitator volume:		220 ml	635 ml	1475 ml
Dimensions (length x diameter):		168 x 71 mm	234 x 99 mm	295 x 127 mm
Weight, incl. reservoir, O ₂ tube and mask:		140 g	215 g	314 g

The maximum frequency of ventilation depends on the insufflation volume used. The rate of expansion of the resuscitator is more than adequate for all frequencies recommended for use.

	Infant:	Paediatric:	Adult:
Pressure limiting valve/ Maximum attainable delivery pressure*:	4.0 kPa (40 cm H ₂ O)	4.0 kPa (40 cm H ₂ O)	4,0 kPa** (40 cm H ₂ O)
Dead space:	< 6 ml	< 6 ml	< 6 ml
Inspiratory resistance without O ₂ :	max. 0.10 kPa (1.0 cm H ₂ O) at 50 l/min	max. 0.50 kPa (5.0 cm H ₂ O) at 50 l/min	max. 0.50 kPa (5.0 cm H ₂ O) at 50 l/min
Expiratory resistance:	0.2 kPa (2.0 cm H ₂ O) at 50 l/min	0.27 kPa (2.7 cm H ₂ O) at 50 l/min	0.27 kPa (2.7 cm H ₂ O) at 50 l/min
Reservoir volume:	300 ml (bag) 100 ml (tube)	2600 ml	2600 ml

*Higher delivery pressure can be obtained by overriding the pressure limiting valve

**Also available without pressure limiting valve


	All versions:
Patient connector:	Outside: 22 mm male (ANSI/ISO) Inside: 15 mm female (ANSI/ISO)
Expiratory connector (for PEEP valve attachment):	30 mm male (ISO)
Forward and backward leak:	not measurable
MediPort:	Standard Luer LS 6
Recommended operating temperature:	-18°C to +50°C
Storage:	Tested at -40°C and +60°C according to EN ISO 10651-4:2002 and ASTM F 920-93

4. Principles of operation

4.1 Resuscitator function


The figures show how the ventilation gas mixture flows into the bag and to and from the patient during manual operation of the resuscitator. The gas flow is similar when the patient is breathing spontaneously through the device.

Adult and Paediatric resuscitator


*Adult version: Also available without pressure limiting valve and manometer port

Infant with bag reservoir

Infant with tube reservoir

4.2 MediPort function

5. Supplementary oxygen

The O₂ reservoir system provides optimal oxygen economy. During compression of the resuscitator oxygen flows through the O₂ inlet tube and further back into the reservoir.

When the resuscitator is released, oxygen is drawn from the reservoir through the inlet valve into the resuscitator. O₂-inlet dimensions is according to DS/EN 13544-2.

Note: If high ventilation pressures are used, higher O₂ flow settings are needed because part of the stroke volume is vented from the pressure-limiting valve. (Adult version: This only applies to the version with pressure limiting valve).

Examples of maximum O₂ percentages, which can be obtained with different volumes and frequencies:

5.1 Oxygen administration for Adult version

O ₂ flow	Ventilation volume (ml) x frequency, I:E ratio = 1:2 Oxygen concentration in %			
Litres/min	250 x 12	600 x 12	750 x 12	1000 x 12
2	74	43	38	34
5	100	76	65	54
10	100	100	100	87
15	100	100	100	100

5.2 Oxygen administration for Paediatric version

O ₂ flow	Ventilation volume (ml) x frequency, I:E ratio = 1:2 Oxygen concentration in %			
Litres/min	40 x 40	100 x 20	200 x 20	400 x 15
1	70	60	40	34
2	100	100	60	47
4	100	100	100	73
6	100	100	100	100

5.3 Oxygen administration for Infant version

O ₂ flow	Ventilation volume (ml) x frequency, I:E ratio = 1:1 Oxygen concentration in %					
	40 x 40		100 x 20		150 x 20	
Litres/ min	Reservoir		Reservoir		Reservoir	
	Bag	10" Tube	Bag	10" Tube	Bag	10" Tube
1	70	70	60	60	47	47
2	100	100	100	100	73	73
4	100	100	100	100	100	100
6	100	100	100	100	100	100
8	100	100	100	100	100	100

Note: For Infant version, use of supplementary oxygen without reservoir attached will limit the oxygen concentration to 60-80% at 15 litres of O₂ /min.

Caution bag reservoir system: The O₂ reservoir bag on the Adult and Paediatric versions are permanently attached to the inlet valve assembly. Do not attempt to disassemble. Tearing may occur. For the Infant version, do not attempt to disassemble the reservoir bag attachment by pulling the bag. Tearing may occur.

6. Operating instructions

6.1 Resuscitation function

Clear mouth and airways using recommended techniques.

Use recommended techniques to position the patient correctly to open the airways and to hold the mask firmly against the face.

It is very important to thoroughly train correct application of the face mask before any attempt is made to use a resuscitator.

Slide your hand (Adult version) or ring and middle finger (Paediatric version) under the support strap. The Infant version does not have a support strap.

Ventilate the patient. During insufflation observe the rise of the patient's chest.

Release the bag abruptly, listen for the expiratory flow from the patient valve and observe lowering of the chest. If continued resistance to insufflation is encountered, check for airway obstruction or correct the backward tilt of the head.

Ventilation without using the support strap (Adult and Paediatric version) can be achieved by turning the bag.

The correct ventilation frequency may vary. Please follow the current ventilation frequency recommended by national or international guidelines.

A pressure gauge can be connected to the manometer port on the top of the patient valve for monitoring ventilation pressure (Adult version: This only applies to the version with manometer port). Remove the cap and connect pressure manometer or the tube of the pressure gauge.

The resuscitator is equipped with an adaptor for the manometer port enabling different types of pressure gauge tubing (Adult version: This applies only to the version with manometer port).


Warning: Use only for monitoring pressure. The cap must always be put on the connector when the pressure is not being monitored.

The pressure-limiting valve is set to open at 40 cm H₂O. If higher pressures are required, place the index finger on the red button while squeezing the bag. The override clip can be used as an alternative to placing the index finger on the red button (adult version: This only applies to the version with pressure limiting valve).

Warning: Never override the pressure-limiting valve unless medical and professional assessment indicate the necessity. High ventilation pressures may cause lung rupture to certain patients. If the pressure-limiting valve is overridden, a manometer must be used to monitor ventilatory pressures to avoid the possibility of lung rupture.

Here the pressure-limiting valve is set to open at 40 cm H₂O.

Here the pressure-limiting valve is closed by means of the override clip.

Remove the outlet cap to fit the Ambu PEEP valve (if required).

Administer oxygen according to medical indications. Refer to "Supplementary Oxygen" (paragraph 5).

6.2 MediPort instruction

Note: AHA 2000 guidelines for administering optimal medication dosages to neonatals by the endotracheal route are not established due to drug absorption variations in patient airways.

Carefully observe patient response to the administered medication.

Administration of volumes of 1 ml fluid or above through the Mediport is comparable with administration directly into an endotracheal tube.

The MediPort has been tested with epinephrine, lidocaine and atropine.

Caution

An increase in the variation of the dosage of medication actually delivered must be expected when administering volumes below 1 ml fluid and without subsequent flushing with appropriate fluid.

Consult your medical director for proper dosing guidelines.

Change to injection directly in the tube if unusually high flow resistance is felt through the MediPort.

Warning

Use only MediPort for drug administration. The MediPort cap must always be put on the MediPort when not administering drugs. Do not attach oxygen supply tubing or any other tubing to this port.

To ensure proper delivery of the entire dosage, the medication port must be flushed after each use.

Syringe with Luer cone

Remove the MediPort cap. Mount the syringe in the MediPort and lock it by turning it 1/4 turn clockwise. Inject drug into the MediPort. Ventilate 5-10 times quickly in succession.

Remove the empty syringe, and replace the MediPort cap.


Syringe with needle

Insert the needle into the middle of the MediPort cap. Inject drug into the MediPort. Ventilate 5-10 times quickly in succession.

Remove the empty syringe.

7. Optional accessory

Ambu PEEP valve

For further information please refer to the directions for use of the Ambu PEEP valve.

Ambu Disposable Pressure Manometer

For further information please refer to the directions for use of the Ambu Disposable Pressure Manometer.

8. Test for correct function

8.1 Resuscitator

Connect a 1.5-2 litre (Adult) or 0.5 litre (Paediatric and Infant) test lung to the patient connector.

Squeeze and release the resuscitator several times and check that the test lung fills. During continuous ventilation, expansion and relaxation of the test bag must be visible. If not, check the inlet valve shutter and the patient valve shutter.

Close the pressure-limiting valve with the override clip (Adult version: This only applies to the version with pressure limiting valve) and the patient connector with a thumb while compressing the resuscitator bag firmly to check tightness and proper valve fitting. (Use a finger to close the pressure limiting valve if the override clip is not present).

Open the pressure-limiting valve by opening the override clip or by removing the finger and repeating the procedure. The pressure limiting valve should now be activated and it should be possible to hear the expiratory flow from the valve.

8.2. Oxygen reservoir system

Resuscitator with reservoir bag attachment

Supply a gas flow of 3 l/min to the oxygen tube. Check that the reservoir fills.

If not, check the integrity of the two valve shutters, or for a torn reservoir or a blocked oxygen tube.

Resuscitator with tube reservoir (Infant version only)

Supply a gas flow of 10 l/min to the oxygen tube. Check that the oxygen flows out at the end of the reservoir tube.

If not, check for a blocked oxygen tube.

8.3. Test for unblocked MediPort passage

Remove the MediPort cap and block the patient connector. Squeeze the bag and listen for the sound of air being pressed out through the MediPort.

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Ambu Ltd.
Ambu B.V.