EN AMBULATOIRE

Ambu® aScope™ 4 RhinoLaryngo Intervention

Ambu

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First experiences with a new flexible single-use rhino-laryngoscope with working channel - a preliminary study

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RÉFÉRENCES

1. Al-Jabir A, Kerwan A, Nicola M, Alsafi Z, Khan M, Sohrabi C, et al. Impact of the Coronavirus (COVID-19) pandemic on surgical practice - Part 1. *Int J Surg.* 2020 ; 79:168-79.
2. Écosse, H.I., 2018. What is the clinical and cost effectiveness of outpatient biopsy for diagnosis of suspicious lesions of the larynx, pharynx and tongue base?. En réponse à une demande du Queen Elizabeth University Hospital, NHS Greater Glasgow and Clyde, [en ligne] Evidence Note (Number 84). Disponible à l'adresse : <http://file:///C:/Users/smpe/Downloads/Evidence-Note-84-Outpatient-Biopsies.pdf> ; [Consulté le 9 avril 2021].
3. Wellenstein DJ, Honings J, Schutte HW, Herruer JM, van den Hoogen FJA, Marres HAM, et al. Cost analysis of office-based transnasal esophagoscopy. *Eur Arch Otorhinolaryngol.* 2019;276(5):1457-63.
4. Schimberg AS, Wellenstein DJ, van den Broek EM, Honings J, van den Hoogen FJA, Marres HAM, et al. Office-based vs. operating room-performed laryngopharyngeal surgery: a review of cost differences. *Eur Arch Otorhinolaryngol.* 2019; 276(11): 2963-73.
5. Wellenstein DJ, Schutte HW, Takes RP, Honings J, Marres HAM, Burns JA, et al. Office-Based Procedures for the Diagnosis and Treatment of Laryngeal Pathology. *J Voice.* 2018;32(4):502-13.
6. Wellenstein DJ, de Witt JK, Schutte HW, Honings J, van den Hoogen FJA, Marres HAM, et al. Safety of flexible endoscopic biopsy of the pharynx and larynx under topical anesthesia. *Eur Arch Otorhinolaryngol.* 2017;274(9):3471-6.
7. Cooper, F., Baig, H., Asimakopoulos, P., Kim W Ah-See, & Shakeel, M. Successful adaptation during COVID-19 pandemic: suspected head and neck cancer biopsies under local anaesthetic. *MedNEXT Journal of Medical and Health Sciences,* 2020;1(1), 15-18.
8. Cancer Research UK. Statistiques sur le cancer de la tête et du cou. Disponible à l'adresse : <https://www.cancerresearchuk.org/health-professional/cancer-statistics/statistics-by-cancer-type/head-and-neck-cancers>. Consulté le 12/02/21.
9. Araghi AS, Harris Y, Kyzas P. An audit to analyse the two-week wait pathway at an oral cancer specialist district general hospital. *Br Dent J.* 2020;228(9):693-699.
10. National Cancer Waiting Times Monitoring Dataset Guidance - Version 11.0. Disponible à l'adresse : <https://www.england.nhs.uk/statistics/wp-content/uploads/sites/2/2020/09/national-cancer-waiting-times-monitoring-dataset-guidance-v11-sep2020.pdf>. Consulté le 25/02/21.
11. Hanna TP, King WD, Thibodeau S, Jalink M, Paulin GA, Harvey-Jones E, et al. Mortality due to cancer treatment delay: systematic review and meta-analysis. *Bmj.* 2020;371:m4087.
12. Sud A, Torr B, Jones ME, Broggio J, Scott S, Loveday C, et al. Effect of delays in the 2-week-wait cancer referral pathway during the COVID-19 pandemic on cancer survival in the UK: a modelling study. *Lancet Oncol.* 2020;21(8):1035-44.
13. Lippert D, Hoffman MR, Dang P, McCulloch TM, Hartig GK, Dailey SH. In-office biopsy of upper airway lesions: safety, tolerance, and effect on time to treatment. *Laryngoscope.* 2015;125(4):919-23.
14. Utility of RhinoLaryngo Intervention scopes in the outpatient setting. Disponible à l'adresse : <https://www.ambu.co.uk/education/educational-videos>. Consulté le 25/02/21.
15. Becker S, Hagemann J, O'Brian C, Weber V, Döge J, Helling K, et al. First experiences with a new flexible single-use rhino-laryngoscope with working channel - a preliminary study. *Laryngorhinotologie.* 2019;98(S 02):11333.

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Ambu A/S