

Release note - Manikin Management Module

V1.1.0

The following has been updated for improved functionality

- Graphs in debriefing did not display data after approximately 2.5 minute.
- When new session was started directly after another, the event icons from the previous session was still displayed on the graph.
- Event icons upload improved.
- Sound files upload improved.
- Added at 3 person CPR training mode.
- The ECC rate shown in a live session is now based on data from the last 2 compressions.
- Improved evaluation on the training results to be more detailed, showing total/percentage values of the performed actions.
- Added "Evaluation" options in Debriefing and print configuration.
- A measure tool is implemented in the graphs and allows to measure seconds between two points by click and dragging.
- Added possibility to ignore initial ventilation. An "Ignore initial ventilation" field is added to all existing algorithms with the default value "No". If initial ventilation is chosen to be ignored, then all ventilations before the first compression is ignored when calculating the ventilation/compression ratio.
- Added a back button on a views in Viewing mode, except the Monitor View.
- Playback of sounds are synchronized to the value of the breathing parameter. All existing default breathing sounds has been updated to contain only one inhale-exhale and will replace the old sounds then the manikin is updated.
- A "Default breathing sound" field is added to all existing ECG signals with the default value depending on the ECG signals breathing parameter.
- Ability to select a sound to play only once added to the sound menu.
- The ECG signal output was changed to lead II instead of lead I.
- The names "master" and "slave" was replaced with "host" and "client".
- The wording "RR" was replaced with "BP".
- All ECG's have been reworked in high resolution. New ECG's added.
- Added Korean language to the system.

Note:

Always save all locally created data and session from the manikin, before performing the update.

With installation of Version V1.1.0 it will not be possible to set back the manikins to a previous version.

Ballerup, April 2018

Earlier releases

V1.0.5

The following reported bug have been updated for improved functionality

Scenarios

- System crash after compression trigger
- Time trigger not counting correct

System

- Configuration page accepts illegal WEP key length
- Invalid network configuration seems to inhibit the network reset button
- Missing text in language files
- Power LED flashes continuously due failure to complete bootup procedure

Control mode

- Defib symbol not shown in ECG signal field
- ECC rate not correct in debriefing table
- ECC and ventilation rate not shown on manikins
- Visibly of Database entries
- Compression / ventilation levels continuously incrementing
- Jump in ventilation graphs

Note:

Always save all locally created data and session from the manikin, before performing the update.

Ballerup, June 2017