








Ambu® Wireless

Next Generation



NEW FEATURES UPDATE V1.1.0

-  Added new and reworked ECG rhythms in higher resolution, compatible with all common defibrillators & patient monitors
-  New possibilities for debriefing added: Individual time measurement in the graphs by clicking and dragging (e.g. hands-off-time, intubation time)
-  Besides the 1, 2 and 4 person CPR – NEW 3 person CPR (training mode) added
-  More detailed evaluation of the training results (view can be set individually in the configuration)
-  Sounds are synchronized to the breathing parameter
-  Ability to play a sound only once in a session
-  Option of initial ventilation

Note: Always save all locally created data and session from the manikin, before performing the update.
With installation of Version V1.1.0 it will not be possible to set back the manikins to a previous version.



OPTIMISATION

- ✓ Improved upload of event icons and sound files
- ✓ Real time ECC rate adapts now faster
- ✓ Workflow optimization: Added a back button on viewing modes (except: patient monitor)
- ✓ ECG signal output changed to lead II instead of lead I
- ✓ Missing graphical data in debriefing were added
- ✓ „Master/Slave“ renamed to „Host/Client“

