



b

Overview of Ambu's Training Devices in relation to content of Latex and REACH Regulated SVHC substances.

According to the European Chemical regulation (REACH) Ambu is obligated to inform about the products containing SVHC substances (Substances of Very High Concern) in the concentrations above 0.1% (w/w).

The table below provides an overview on all Ambu Training Devices in relation to content of REACH Regulated SVHC substances in concentrations above 0.1%. The table also includes information on which products contain Natural Rubber Latex.

Product	Contents of SVHC above 0.1% (w/w) and description of substance	Natural Rubber Latex	Symbol on label
Emergency Care and Training			
BLS:			
Ambu® Man	Yes, DEHP	No	
Ambu® Baby Manikin	Yes, DEHP	No	
Ambu® Defib Trainer System	Yes, DEHP	No	
Ambu® CPR Pal	Yes, DEHP	No	
Ambu® Man Compression	No	No	
Ambu® Man Basic	Yes, DEHP	No	
Ambu® SAM	Yes, DEHP	No	
Ambu® Junior	No	No	
Ambu® Man Wireless - Next Generation	Yes, DEHP	No	
Ambu® Man Instrument	Yes, DEHP	No	
Ambu® Man Defib – Next Generation	Yes, DEHP	No	
Ambu® Man School	No	No	
Ambu® Patented Hygienic System	Yes, DEHP	No	
Ambu® CPR software	No	No	
ALS:			
Ambu® Airway Man	Yes, DEHP	No	
Ambu® Airway Management Trainer	No	No	
Ambu® IV Trainer	Yes, DEHP	Yes	
Ambu® Kit i.o. Access	No	No	
Ambu® Man Airway Wireless – Next Generation	Yes, DEHP	No	
Ambu® Man Airway Instrument	Yes, DEHP	No	
Ambu® Man Advanced IV – Next Generation	Yes, DEHP	Yes	
Ambu® Man Advanced – Next Generation	Yes, DEHP	No	
Ambu® Cardiac Care Trainer	Yes, DEHP	Yes	

Ambu A/S
Baltorpbakken 13
DK-2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu@ambu.com
www.ambu.com
CVR. nr. 63644919