

# Ambu® AuraGain™



- Le nouveau niveau de sécurité et d'efficacité



—  
ana d'acc s astri ue  
ossibi ité d'intubation  
ourbe anatomiT□

□D□inflatable cuff  
o D□Phthalate

Nouvelles  
tailles  
Pédiatrique





## Specifications

| Ambu AuraGain Pediatric |       |         |          |          |
|-------------------------|-------|---------|----------|----------|
| Mask size               | #1    | #1½     | #2       | #2½      |
| Packaging color code    | ●     | ●       | ●        | ●        |
| Patient weight          | <5 kg | 5-10 kg | 10-20 kg | 20-30 kg |
| Max. ETT size           | 3,5   | 4,0     | 5,0      | 5,5      |
| Max gastric tube size   | 6 Fr  | 8 Fr    | 10 Fr    | 10 Fr    |

Visit our website at: [ambu.com/AuraGain](http://ambu.com/AuraGain)

**Ambu A/S** • Baltorpbakken 13 • DK - 2750 Ballerup • Denmark

Tel. +45 7225 2000 • Fax +45 7225 2053